

KI0071 HPLC Analysis Report

Sample: BIBW2992

Sample ID: KI0071

Buffer A: 0.1% TFA in 100% water (v/v)

Buffer B: 0.1% TFA in 80% acetonitrile + 20% water (v/v)


Gradient: 10–100% Buffer B in 30min

Flow: 1ml/min

Wavelength: 220nm

Column: Venusil C18 5um 4.6*250mm

Chromatogram



Signal	Peak#	RT(min)	Area%	Area(mAU*s)	Height (mAU)
	1	15.300	0.2977	50526	8208
	2	15.447	99.05	16810094	1504354
	3	15.983	0.6551	111183	19816
Total			100	16971803	1532378