

HPLC Analysis Report

Sample: Dutasteride

CAS No.: 164656-23-9

Batch No.: KB20140604

Mfg. Date: June 04, 2014

Method Info: 30%-90% Buffer B 0-25 MIN

Buffer A: 0.04% TFA in 100% water(v/v)


Buffer B: 0.04% TFA in 100% acetonitrile (v/v)

Column: Alltima C18,4.6*250mm,5um

Flow Rate: 1mL/min

Wavelength: 220 nm

Date: June 10, 2014



Signal Peak#	RT(min)	Width(min)	Area(mAU*s)	Height(mAU)	Area%
1	4.139	0.1136	9.03484	1.16712	0.1901
2	18.622	0.1613	4743.84766	453.63712	99.8099
Total			4752.88250	454.80423	100.0000

Analysis Result:

Purity: 99.81%