

KI0307 HPLC Analysis Report

Sample: Pazopanib

Sample ID: KI0307

Buffer A: 0.1% TFA in 100% water (v/v)

Buffer B: 0.1% TFA in 80% acetonitrile + 20% water (v/v)


Gradient: 10-100% Buffer B in 30min

Flow: 1ml/min

Wavelength: 220nm

Column: Venusil C18 5um 4.6\*250mm

Chromatogram



Signal	Peak#	RT (min)	Area%	Area (mAU*s)	Height (mAU)
	1	15.960	0.1708	43253	8663
	2	16.539	99.83	25287443	4421663
Total			100	25330696	4430326