

LCMS Analysis Report

Sample: Raloxifene(Evista), CAS: 84449-90-1

Buffer A: 0.04% TFA in 100% water(v/v)

Buffer B: 0.04% TFA in 100% acetonitrile (v/v)


Gradient: 10-80% Buffer B in 25 min


Flow: 1ml/min Wavelength:220nm

Column: Alltima C18,4.6*250mm,5um

Date: May. 27, 2014

Chromatogram





Peak#	MSBase(+)	Height(mAU)	Area(%)
1	ND	5.88	0.526
2	474.1	730.24	99.474